

# HiSeq™ 2000 Sequencing System

Redefining the trajectory of sequencing.

## HiSeq 2000 Highlights

- High Accuracy and Unprecedented Output:**  
Generate up to 200 Gb per run with the highest yield of data greater than Q30, 2 × 100 bp read length, up to 25 Gb per day, two billion paired-end reads per run.
- Breakthrough User Experience:**  
Easily set up runs with simplified library prep, automated clonal amplification, pre-configured, plug-and-play reagents, simple flow cell loading, touch screen-enabled user interface, and integrated paired-end fluidics.
- Unmatched Cost-Effectiveness:**  
Unrivalled output and ease of use provide the lowest overall operating cost. Sequence two human genomes at greater than 30× coverage in 8 days. Process 200 gene expression samples in a single run at a per sample price less than microarrays.

## Sequence at a Scale Never Before Possible

The HiSeq 2000 sequencing system (Figure 1) offers unprecedented output and a breakthrough user experience. Leveraging Illumina's proven and widely-adopted, reversible terminator-based sequencing by synthesis chemistry in combination with innovative engineering, HiSeq 2000 delivers the industry's highest sequencing output and fastest data generation rate. Human interaction design features and the easiest sequencing workflow set a new standard for simplicity and user experience. With unmatched cost-effectiveness, HiSeq 2000 is the first commercially available sequencer to enable researchers to obtain ~30× coverage of two human genomes in a single run for under \$10,000 (USD)\* per sample.

## Unprecedented Output

HiSeq 2000 makes it possible for individual labs to take on the largest and most complex sequencing studies at the lowest cost. With cutting-edge scanning and imaging technology, clusters on both surfaces of the flow cell can be sequenced, dramatically increasing the number of reads, sequence output, and data generation rate. The ultra-high output and speed of HiSeq 2000 now make it possible to sequence two human genomes at ~30× coverage simultaneously, up to 200 gene expression samples in a single run, or one billion tags per day for short-read applications.

\*Reagent costs at list price.

## Breakthrough User Experience

Innovative design features make HiSeq 2000 the easiest-to-use next-generation sequencing system (Figure 2). Flow cells are loaded on the vacuum-controlled loading dock. Pre-configured, plug-and-play reagents sufficient for up to 200 cycles drop into racks in the machine's chiller compartment, requiring only two minutes of hands-on time. A simple touch screen user interface, including on-screen, step-by-step instructions with embedded multimedia help, simplifies run setup. Real-time progress indicators provide at-a-glance status, and remote monitoring allows a single user to check progress on multiple systems from any browser or internet-enabled phone. HiSeq 2000 can be operated in single or dual flow cell mode, offering unmatched experimental flexibility and instrument scalability. Independently-operable flow cells allow applications requiring different read lengths to run simultaneously.

Figure 1: HiSeq 2000



Illumina's HiSeq 2000 sequencing system enables individual labs to take on larger and more complex studies, including routine human genome sequencing.





## Accelerate Your Research with HiSeq 2000

HiSeq 2000 redefines the trajectory of sequencing by combining innovative engineering with proven sequencing by synthesis chemistry to set new standards for output, simplicity, and cost-effectiveness. With HiSeq 2000, the ability to process larger numbers of samples and to decode larger and more complex genomes means that virtually any sequencing project is now within reach.

## Learn More

For more information about HiSeq 2000 and Illumina sequencing, visit [www.illumina.com/hiseqnow](http://www.illumina.com/hiseqnow).

## HiSeq 2000 Information

### HiSeq 2000 Performance Parameters\*

Read Length	Run Time	Output
1 × 35 bp	~1.5 days	26–35 Gb
2 × 50 bp	~4 days	75–100 Gb
2 × 100 bp	~8 days	150–200 Gb

\*Sequencing output generated with a PhiX library and cluster densities between 260,000–347,000 clusters/mm<sup>2</sup> that pass filtering on a HiSeq 2000.

### Throughput

Up to 25 Gb per day for a 2 × 100 bp run.

### Reads

Up to one billion clusters passing filter, and up to two billion paired-end reads.

### Performance

HiSeq 2000 provides the greatest yield of perfect reads and bases greater than Q30:

Greater than 90% bases higher than Q30 at 2 × 50 bp\*\*

Greater than 85% bases higher than Q30 at 2 × 100 bp\*\*

\*\*Human genome at supported cluster densities

### Service and Support

Illumina will ensure that your HiSeq 2000 is properly installed and qualified, and will provide ongoing maintenance and service. This industry-leading support is available in North America, Europe, and Asia.

## HiSeq 2000 System Specifications With Monitor and PC

### Catalog No. SY-401-1001

#### Instrument Configuration

CE marked and ETL listed instrument (HiSeq 2000 only)  
Computer and touch screen display  
Installation setup and accessories  
Data collection and analysis software

#### Instrument Control Computer

Base Unit: 2x Intel Xeon X5560 2.8 GHz CPU  
Memory: 48 GB RAM  
Hard Drive: 4x 1.0 TB 7200 RPM SATA  
Operating System: Windows Vista  
Note: Computer specifications will be regularly upgraded.  
Contact your local account manager for current configuration.

#### Operating Environment

Temperature: 22°C ± 3°C  
Humidity: Non-condensing 20%–80%  
Altitude: Less than 2,000 m (6,500 ft)  
Air Quality: Pollution degree rating of II  
Ventilation: Maximum of 4,000 BTU/h  
For Indoor Use Only

#### Laser

532 nm, 660 nm, 700–840 nm, 650 nm (barcode reader)

#### Dimensions

WxDxH: 118.6 cm × 76.0 cm × 94.0 cm  
(46.7 in × 30.0 in × 37.0 in)  
Weight: 221.4 kg (488 lbs)  
Crated Weight: 312 kg (688 lbs)

#### Power Requirements

100–240V AC 50/60Hz, 20A, 1500W

Illumina provides a region-specific uninterruptible power supply with the HiSeq 2000 instrument.

## HiSeq 2000 Kits and Accessories

Application	Catalog No.
HiSeq 2000 Sequencing System	SY-401-1001
cBot Clonal Amplification System	SY-301-2002
Paired-End DNA Sample Prep Kit	PE-102-1001
Genomic DNA Sample Prep Kit	FC-102-1001
HiSeq Paired-End Cluster Generation Kit	PE-401-1001
HiSeq Single-Read Cluster Generation Kit	GD-401-1001
HiSeq Sequencing Kit (200 cycles)	FC-401-1001
HiSeq Sequencing Kit (50 cycles)	FC-401-1002

### FOR RESEARCH USE ONLY